



Quality control for scRNAseq data

Orr Ashenberg, Jacques Serizay, Fabrício Almeida-Silva

November, 2024




Lecture topics

- Interacting with Seurat objects
- Quality control, normalization, and feature selection starting from raw count or expression matrices
- Next step will be dimensionality reduction, clustering, and visualization

Determining cell type, state, and function

Quality control



Normalization

Feature selection

Dimensional reduction

Cell-cell distances

Unsupervised clustering

Interacting with Seurat objects

- Seurat object is used for to store 10x data and perform analysis
 - Count matrices for different assays are stored (gene expression, protein expression, chromatin accessibility, etc...)
 - Counts are stored as: counts (raw), data (normalized), scaled data (centered and scaled) in sparse matrices when possible
 - Metadata describes individual cells and genes
 - Functions for analysis (quality control, normalization, feature selection, dimensional reduction, cell-cell distances, unsupervised clustering)

<https://github.com/satijalab/seurat/wiki>

https://satijalab.org/seurat/essential_commands.html

Interacting with Seurat objects

```
> gdata
```

```
An object of class Seurat
35633 features across 2000 samples within 2 assays
Active assay: RNA (33633 features)
 1 other assay present: integrated
 2 dimensional reductions calculated: pca, umap
```

Seurat object

```
> gdata[['RNA']]@data[1:5,1:5]
```

```
5 x 5 sparse Matrix of class "dgCMatix"
```

	D2ex_5	D2ex_6	D2ex_7	D2ex_11	D2ex_13
A1BG-AS1
A1BG
A1CF
A2M-AS1
A2ML1	.	.	.	1.226772	.

Accessing count slot from RNA assay

```
> gdata[[1]][1:5, 1:5]
```

	orig.ident	nCount_RNA	nFeature_RNA	tech	integrated_snn_res.1
D2ex_5	D2ex	5745.867	2548	celseq	4
D2ex_6	D2ex	6883.692	2619	celseq	6
D2ex_7	D2ex	7460.202	3043	celseq	5
D2ex_11	D2ex	8330.644	3465	celseq	5
D2ex_13	D2ex	3891.960	1962	celseq	6

Accessing cell metadata

```
> gdata <- ScaleData(gdata)
```

```
Centering and scaling data matrix
```

Running analysis function

```
|=====| 100%
```

Loading data into a Seurat object

```
gcdata <- CreateSeuratObject(counts = celseq.data)
```



counts matrix

	Cell 1	Cell 2	Cell 3	...	Cell 5K
Gene 1	3	0	1		2
Gene 2	0	2	0		1
Gene 3	1	0	3		5
...					
Gene K	14	7	1		0
...					
Gene 25K	0	13	1		0

Storing counts data in dense vs sparse format

2D Arrays

cells								genes
1	2	3	4	5	6	7	8	
0	0	0	0	0	0	0	0	
0	1	0	0	0	0	0	0	
0	0	0	0	0	2	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	0	
0	0	0	0	0	0	0	3	6

Dense matrices


Coordinate List

2	2	1
6	3	2
8	6	3

Sparse matrices


There are many quality control filters for genes and cells

Complexity =
Number of genes
detected in a cell



There are many quality control filters for genes and cells

- We filter cells based on technical or biological parameters.




Filtering with combinations of quality control filters

Low nUMI and high % mitochondrial-
Cells captured but lost a lot of the
mRNA, and the mitochondrial genes
were protected and retained

**High nUMI & low nGene ratio – low
quality library or capture rate**

High nUMI & high nGene – doublets




Appropriate quality control filters vary with platform and cell types

- Different platforms set different expectations
 - e.g. Smart-Seq2 often yields more genes detected per cell than 10x Chromium.
- Different cell types set different expectations
 - Immune cells normally have fewer genes detected per cell than non-immune cells
 - Malignant cells normally have more genes detected per cell than non-malignant cells

A classifier for low-quality cells

C



“The pipeline takes advantage of a highly-curated set of generic features that are incorporated into a machine learning algorithm to identify low quality cells.”

A classifier for low-quality cells

SVM feature weights






What are caveats
to this approach?

Other quality control filters for genes and cells

- Doublets

- number of genes 
- number of UMIs 
- percentage of mitochondrial gene expression 


- Ambient RNA and empty droplets

- number of genes 
- number of UMIs 
- percentage of mitochondrial gene expression 

- Barcode swapping

Cell doublets can be misleading

Because of the setup, it is possible that two or more cells can enter the same droplet. Studies estimate doublet frequency through a “mixed-species” experiment



The doublet frequency is positively correlated with throughput


Detecting cell doublets with Scrublet

Scrublet (Single-Cell Remover of Doublets)

Singlets



Detecting cell doublets with Scrublet




Detecting empty drops containing ambient RNA – manual

Look for transcripts expressed in unexpected cell types and remove those genes from all subsequent analysis


- e.g. hemoglobin expressed in a T cell

Detecting empty drops containing ambient RNA – automatic

EmptyDrops (distinguish cells from empty droplets)



Detecting barcode swapping in multiplexed samples




“Barcode swapping is a phenomenon that occurs upon multiplexing samples on the Illumina 4000 sequencer. Molecules from one sample are incorrectly labelled with *sample* barcodes from another sample, resulting in their misassignment upon demultiplexing.”

”Specifically, we considered molecules across **multiplexed samples that contain the same combination of unique molecular identifier, cell barcode, and aligned gene.**”


Tool to detect empty drops, correct ambient RNA and barcode swapping

CellBender remove ambient background and barcode swapping via deep learning


Detecting empty droplets



Correcting ambient RNA and barcode swapping




Determining cell type, state, and function





Normalizing gene expression in each cell

- Why normalize gene expression within a cell?
 - cells are sequenced to different depths (technical)
 - cells of different type have different amounts of mRNA (biological)
 - there are typically extreme values in distribution of gene expression
 - more highly expressed genes are more variable




Normalizing gene expression in each cell

- How to normalize
 - Gene expression measurements for each cell are normalized by the total gene expression or median gene expression
 - Gene expression values then scaled to sum to 10,000 (typically), and then log-transformed.




Is standard normalization appropriate?

Reassessing the idea that droplet scRNA-Seq is zero-inflated




- “Droplet scRNA-seq is not zero-inflated.” Svensson, Nature Biotechnology (2020)
- “Feature selection and dimension reduction for single-cell RNA-Seq based on a multinomial model.” Townes et al. Genome Biology (2019)
- “Normalization and variance stabilization of single-cell RNA-seq data using regularized negative binomial regression.” Hafemeister et al. Genome Biology (2019)


Is standard normalization appropriate?



(a) UMI counts




(b) counts per million (CPM)



(c) log of CPM

Example of how current approaches to normalization and transformation artificially distort differences between zero and nonzero counts. **a** UMI count distribution for gene ENSG00000114391 in the monocytes biological replicates negative control dataset. **b** Counts per million (CPM) distribution for the exact same count data. **c** Distribution of $\log_2(1+CPM)$ values for the exact same count data


Determining cell type, state, and function



Identify highly variable genes

Cells are in ~20,000 dimensional space (one dimension for each gene)

- many genes are lowly detected or noisy measurements




- variable genes contain the biological signal we are interested in

Calculating gene signatures

Relying on capturing a specific gene is not robust, but relying on a set of genes (signature) is much more stable!



Gene signature example: cell cycle markers



variability of individual genes

